

A COMPREHENSIVE GUIDE TO BISULFITE-CONVERTED DNA AMPLIFICATION

Over the years, bisulfite conversion has become the most widely used method for DNA methylation analysis. It is the most convenient and effective way to map DNA methylation to individual bases. As the first step in numerous downstream analysis techniques, it is absolutely critical that the bisulfite conversion process is understood and performed well.


Our scientists at Zymo Research have spent extraordinary amounts of time studying and optimizing bisulfite conversion as well as supporting our customers. Today we see fewer questions arise from the bisulfite conversion process, but more questions related to analyzing bisulfite converted DNA in the diverse types of downstream applications.

In this guide, we've pulled together key points and considerations that should aid researchers of any experience level to better understand key factors that can impact downstream analysis.

TIPS FOR SUCCEEDING DOWNSTREAM OF BISULFITE CONVERSION

Since DNA methylation analyses are just getting started after bisulfite conversion, it's useful to first understand what just happened to your sample during the conversion.

Bisulfite conversion is a relatively harsh process that will dramatically change the chemical and physical properties of your DNA sample. During the process, your genomic DNA transitions from large and stable double-stranded molecules to an assortment of randomly fragmented single-stranded molecules. Additionally, most of the cytosine bases have been converted to uracil.



Basically, your DNA does not resemble DNA anymore and has gone through a major overhaul, so you will have to make some adjustments to your subsequent experimental steps.

CHECKING THE QUALITY OF BISULFITE CONVERTED DNA

Just as with any other experiment, you will want to assess the quality and quantity of your sample before moving forward. A couple of common ways to accomplish this include:

- UV Spectrophotometry
- Agarose Gel Electrophoresis

Keep in mind that it's necessary to make a few adjustments to each of these methods when working with bisulfite converted material.

UV-Quantification of Bisulfite Converted DNA

When using a spectrophotometer, use a value of 40 µg/mL for Ab_{260nm} = 1.0, because the converted DNA now more resembles RNA.

At first, many researchers notice their recovery seems very low. There are two reasons for this:

1. It is low. It is possible to lose sample during the bisulfite treatment process, particularly if degraded DNA is going into the reaction, which can lead to poor recovery.
2. RNA contamination inflated your initial quantification. During the bisulfite conversion and clean-up processes, RNA is removed so your follow up quantification will appear lower.

In either case, do not fear; the recovered material is generally still sufficient, particularly if a PCR step is planned downstream. Starting DNA that is RNA-free and intact will typically yield the best results.


Gel Analysis of Bisulfite Converted DNA

When using an agarose gel to analyze the converted DNA, use a 2% gel with a 100 bp marker. When using gels it can be necessary to load up to 100 ng of sample to visualize the DNA.

Many researchers worry when at first they can't see any band in their gel. This is normal as most of the DNA is single stranded.

To address this, cool the gel for a few minutes in an ice bath. This will drive some base pairing between the single-stranded molecules so the ethidium bromide can intercalate and illuminate your recovered material.

The recovered DNA will usually smear from 1,500 bp down to 100 bp.



The Challenges of Bisulfite-Converted DNA. DNA fragmentation and limited base pairing make it difficult to see the DNA following bisulfite conversion. Samples of bisulfite converted human genomic DNA were run in a 2% agarose gel with a 100 bp marker visualized normally (left panel) and after chilling in an ice bath (right panel)

GETTING STARTED WITH DOWNSTREAM BISULFITE ANALYSIS

There is a tremendous range of bisulfite-based methods out there, but we'll focus on some of the most widely used in labs today (i.e., PCR-based).

PCR-Based Analysis of Bisulfite Converted DNA

Numerous bisulfite-based analysis approaches employ PCR amplification, and some of the more widely used of these methods include:

- Methylation Specific PCR (MSP)
- Bisulfite Sequencing
- Bisulfite Pyrosequencing

The PCR step in these methods is also the most error prone. Again, most of the challenges of amplifying bisulfite converted DNA relate to the fact that the sample's physical traits have changed, including:

- The converted DNA has become highly fragmented
- The strands are no longer complementary
- The molecules have lost most of their cytosine composition, turning it into more of a three-base containing nucleic acid

Although, these conditions do not make for an easy PCR experience, a few practices can mitigate the impact posed from this challenging template including:

- Designing optimized primers
- Using a "hot start" polymerase
- Optimizing annealing conditions

OPTIMIZING BISULFITE PRIMER DESIGN

Optimizing primers for bisulfite amplification can be challenging at first, but if you follow a few guidelines, you should be on your way to successful results relatively quickly. There are two types of primer design approaches depending on your application Bisulfite PCR primers, which will amplify a region regardless of methylation status, and Methylation Specific PCR (MSP) primers which are designed to amplify only when a target locus is either methylated or unmethylated:

Primer Design for Bisulfite PCR

In this scenario, primers are used, as they are in other PCR settings, to amplify the template, most often for analysis by another downstream methodology (e.g. bisulfite sequencing, mass spectrophotometry) that will determine the DNA methylation status.

Featured Tool: Bisulfite Primer Seeker

With over a decade of experience in bisulfite conversion technology, Zymo Research has simplified the tedious process of bisulfite primer design with the [Bisulfite Primer Seeker](#). What sets this program apart from other available options is its ability to design primers from particularly CG-rich templates. The Bisulfite Primer Seeker also stands out by supplying multiple primer options for amplicons that span different regions within your sequence of interest.

Bisulfite Primer Seeker

Sequence

CCCTACAAACCGCTCGTAGAATTCGTGCTCGGCCCTCGTAGTGCGGCCCTCAGTCGCGT
TCCCGAGTAGAGGCGACAGGCGGCGGACACACAGGACAGGGCCCCGTCACCCCTCC
GCAGGCTCCACCAACAGGTATCGCTGGGTCCACTCGGGCCGGAACACAGAGCCTCG
TCGACTTCCATCCTAATCCTCTGGGCGTCATCCACATTCTGCGGAGGCCACAGGAGAA
CGCCGAGCTCTGAGGCCGCGAGCAACGGGGAGGAGCCGACGAACGCCAAGCACCC
CTGGGTTCTGCGCCAGATCACCCGAGCAGGGGTGCACACGCTCCGCCAGGGCACCT
CCGTTCTATTGGCTGGAGATTTCGTATTCCTCAGCTCTCAGAAATAGGCAATAGGAA
GAGCGGACCGTGAACTTTGACGCGCAAGCGCTTGCTTCTAGCCTGGTGTCGGGCCG
CTGTTTAGGGCGGGACAGAGATCCAGGAAGTACGCTGAAAAGCGCTGACTGGCATTC
ATGCTGCC

1000 bp

Parameters

* **Primer Length** From 24 To 38
* **Product Length** From 100 To 350
* **Tm** From 55 To 66
☒ **Allow 1 CpG in the first 1/3 of primer**
* **Email:**
Job Name:
email me my primers

Preview

☒ Line # ☐ BS Conversion ☐ Reverse Complementary

1	GACATGAAGAGCCCCAACTGAAGTGATGAGTCAGGTTGATAATCGTCAGGAAGATCTTCTGCACCTCCAGGTGCGGCCGACGCCACGGATGAAGAC
101	CAGCAGGTAGTTCTCATAGGCCACAAAGCCTGGTCATCCAAAGCAAGGAATAGGCCTTAAAGTCCCTGGCTCGTTAAACAGCTGGCTGCGAGGTTG
201	CTGTCAATGCTCAGAATCCTCTGCGGTGTATCTCTGGAGACAAGTCAATGCCCTCCAGGACGCCGAGTGATCCGGCAGCACCTCTCTCAGCAACACCT
301	CCATGCCTGGTGTACGAAGTCACCTCACCCAGCCGCGACCTTCAAGGCCAAGAAGCGCAGAGGCCGAGGCTGCCGCGTCTCTCTCTCTCCG
401	AGTGAGCACGCGGCCAAAGACATGTCACTCTGCGCAGACGCTCCACAGGGCGCGCTTCTCTATCCCTACAAACCGCTCTAGAAATCTGTCTCG
501	GCCTCTGTAGTGGCCTCACTGCGGTTCCGAGTAGAGGCGACCAAGCGCGACACACAGGACAGGGCCCCGTCACTCCCTCCAGGCTCCACCA
601	GGTATCGCTGGTCCACTCGGGCCGGAACACAGAGCCTCGTCACTTCCATCCTAATCTCTTGGGCGTCATCCACATTCTGCGAGGCCACAGGAGAA
701	CGCCGAGCTCTGAGGCCGAGCAACGGGAGGAGCCGACGAACCAAGCACCCCTGGGTTCTGCGCAGATCACCAGCAGGGGTGACACAGCT
801	CGCCAGGGCACCTCGTTCCTATTGGCTGGAGATTCTGATTCCTCAGCTCTCAGAAATAGGCAATAGGAAGAGCGGACCTGAACCTTTGACCGCA
901	AGCGGCTTGCTTCTAGCCTGGTGTGCGGCTCTGTTAGGGCGGACAGAGATCCAGGAAGTACCTGAAAAGCGCTTGAATGCTGCTGCTGCTG

CpG

Screenshot of Bisulfite Primer Seeker
(www.zymoresearch.com/tools/bisulfite-primer-seeker)

Unlike normal PCR where your template contains 4 bases, bisulfite PCR primers need to be longer to contend with the loss of most of the cytosines. Shoot for primers of at least 26-30 bases of sufficient specificity. Also, since your template will be highly fragmented, we recommend targeting an amplicon size between 150-300 bp.

Next, you will want to make sure that the primers avoid CpG sites, or if that is not possible; try to locate them at the 5' end of the primer with a mixed base at the cytosine position.


When designing the primers, it is a good exercise to 'convert' the DNA first – changing all of the non-CpG cytosines to uracils. It is also important to note that only one strand of the template will be amplified by a given primer set. Because the DNA strands are no longer complementary, only the reverse primer will actually bind to the target DNA. The forward primer will only have a complementary template once generated from the reverse primer. Usually, 35 to 40 cycles are required for successful amplification.

Primer Design for Methylation Specific PCR (MSP)

Unlike primers for bisulfite PCR, MSP primers serve the dual functions of amplifying bisulfite converted DNA and interrogating the DNA methylation status at specific CpG sites.

In addition to observing the above practices for standard bisulfite amplification, MSP primers should be designed taking the following guidelines into account:

1. CpG site(s) of interest must be included in the primers.
2. The CpG(s) should be located at 3' end of the primers.
3. Two primer sets are required for each amplicon that will be interrogated
 - A "Methylated" primer set that contains cytosine (C) at CpG sites.
 - And, a "Non-methylated" primer set that uses thymine (T) in place of C at CpG locations.



Flowchart of primer design for bisulfite PCR and Methylation Specific PCR (MSP). (A) Following bisulfite treatment, the two converted strands of the DNA template are no longer complementary. (B) Primers for Bisulfite PCR are designed for subsequent sequencing and analysis of cytosines within the amplicon. CpG sites within the primers should be avoided or located at the 5'-end with a mixed base at the cytosine position (Y= C/T, R= G/A). Sequencing data is commonly represented by a "lollipop" plot where closed circles represent methylated cytosine positions and open circles non-methylated ones. (C) Primers for Methylation Specific PCR (MSP) are designed to target and assess the methylation status at specific CpG sites. CpG sites within the primers must be located at the 3'-end to increase their specificity to methylated (M) or non-methylated (U) templates. Completely methylated or non-methylated templates will generate a single amplicon from only their representative primer set following MSP. Samples with mixed methylation, will be amplified by both primer sets.

USE A HOT START POLYMERASE

When amplifying bisulfite-converted DNA, it is ideal to use a hot start polymerase. Non-specific amplification is relatively common with bisulfite-treated DNA due to its AT-rich nature.

A hot start polymerase will minimize primer dimerization and the non-specific amplification products that go along with it, better than regular Taq polymerase.

TURN UP THE HEAT WITH HIGHER ANNEALING TEMPERATURES

At Zymo Research, we have found that higher annealing temperatures are key for good amplification, and work outside our lab, including that of M.D. Anderson's Lanlan Shen,

([Biotechniques, 2007](#)) also confirmed that raising the annealing temperatures was instrumental in improving PCR efficiency from bisulfite converted templates.

The good news is that the longer primers suggested for bisulfite primer design, not only help with specificity, but also bring your melting temperature up. It's also best to include as many guanines as possible in the primer regions to increase the melting temperature of the primer to the recommended 55-60 °C range.

When first trying out a new primer set, we suggest running an annealing temperature gradient. This will help you identify the ideal annealing temperature to use in your experiments. Troubleshooting your primers before starting your PCR analysis will save you time and eliminate any non-specific results.

CHOOSING WISELY: WHICH BISULFITE CONVERSION APPROACH IS RIGHT FOR YOU?

A few years ago, there wasn't much to choose from in terms of bisulfite conversion options. Today, there are many options, which can also present some challenges.

Let's run through a few considerations and research scenarios that will impact which type of bisulfite-treatment kit you'll want to use:

- Researchers New to Bisulfite Conversion Technology
- Researchers Seeking Speed & Convenience
- Researchers Working with Cell and Tissue Samples
- Researchers Working with Small Amounts of Precious Samples
- Researchers Seeking Validated Kits for Use in Commercial and Academic Protocols

Spin Column ≥10 µl



96-Well ≥15 µl



MagBead ≥25 µl



Available Kit Formats. All of the bisulfite conversion kits are available in the indicated formats, with the indicated elution volumes. Single columns are ideal for routine microcentrifuge based processing of relatively few samples. 96-well spin plates allow increased throughput for larger studies and sample numbers. Magnetic beads are suited for high-throughput automated handling of large sample sets.

Researchers who are New to Bisulfite Conversion

If you are a new user, and starting a bisulfite conversion project from scratch, you will need a bisulfite conversion kit that's ready to go right out of the box. The EZ DNA Methylation-Lightning™ Kit sports a fast, highly robust and simple protocol, along with pre-mixed reagents, making it a great all-around performer. The ease and efficiency of the Lightning kit is ideal for scientists just entering the realm of DNA methylation.

Research Situation: New Bisulfite Users
Kit Recommendation: EZ DNA Methylation-Lightning™ Kit

Researchers Who Seek Speed and Convenience

Many bisulfite conversion techniques use an incubation time of several hours, or even overnight, and require that reagents be prepared and set up ahead of time. The EZ DNA Methylation-Lightning™ kit is designed for researchers who crave simplicity and don't have time to spare.

The Lightning Kit uses a pre-mixed conversion reagent that is ready-to-use right out of the box, plus it has been optimized to completely convert DNA in about 1 hour.

Research Situation: Seeking Speed and Convenience

Kit Recommendation: EZ DNA Methylation-Lightning™ Kit

Researchers Working with Cell and Tissue Samples

Most bisulfite conversion protocols will generally require you to first isolate the DNA from your sample, and then move ahead with the bisulfite conversion. This is because cellular debris can reduce the efficiency of the conversion reaction. The downside to any purification step though, is that recovery is never 100%, so you can lose precious sample.

To address this, the EZ DNA Methylation-Direct™ Kit combines a unique proteinase digestion step with the conversion reaction enabling robust conversion without having to first perform a DNA purification step. This maximizes the recovery of the DNA template for bisulfite conversion.

Research Situation: Direct Input of Cellular or Tissue Samples

Kit Recommendation: EZ DNA Methylation-Direct™ Kit


Researchers Working with Small Amounts of Precious Samples

Another challenging situation researchers often face is having very little sample to work with (e.g., FFPE, LCM, Needle core biopsies).

Moreover, our scientists have developed a specialized conversion reagent in the EZ DNA Methylation-Direct™ Kit for ultra-low sample inputs, whether the DNA has been purified or not.

Research Situation: Small Sample Inputs

Kit Recommendation: EZ DNA Methylation-Direct™ Kit



Researchers Seeking Validated Kits for Use in Commercial Platforms

Our bisulfite conversion systems are recommended by a number of commercial platforms due to their reliable performance.

Illumina Infinium®

The HumanMethylation450 BeadChip array allows interrogation of DNA methylation status at around 450,000 unique CpGs in a single experiment. R&D scientists at Illumina validated several commercially available kits and found the EZ DNA Methylation™ Kit to be highly compatible with their protocol.

Kit Recommendation: EZ DNA Methylation™ Kit

Sequenom MassARRAY®

MassARRAY® is a MALDI-TOF based system that can measure individual methylation ratios for CpGs within a target sequence to determine relative methylation status of up to several hundred CpGs in multiple samples. The team at Sequenom found that the EZ DNA Methylation™ Kit enabled reliable and consistent results on the MassARRAY® platform.

Kit Recommendation: EZ DNA Methylation™ Kit

Agilent SureSelect Target Enrichment

SureSelect XT™ Methyl-Seq uses sample enrichment techniques combined with bisulfite sequencing to analyze over 3.7 million individual CpG sites for their methylation status at one time. The system targets promoters, CpG islands, and known differentially methylated regions. Scientists at Agilent have shown that EZ DNA Methylation-Gold™ kits have demonstrated very high performance in this system.

Kit Recommendation: EZ DNA Methylation-Gold™ Kit

Epicentre EpiGnome™ Kit

The EpiGnome™ methyl-seq kit creates whole genome bisulfite sequencing libraries from as little as 50 ng of genomic DNA. EpiGnome™ follows a unique “post-bisulfite conversion” library construction method, which yields highly diverse libraries with uniform coverage. Epicentre’s development group recommends the EZ DNA Methylation-Lightning™ kit for use with their product.


Kit Recommendation: EZ DNA Methylation-Lightning™ Kit

Roche NimbleGen SeqCap™ Epi Enrichment System
SeqCap Epi Enrichment enables the targeting of selected genomic regions from bisulfite treated genomic DNA, which can act as a screening tool to identify specific regions in the genome for methylation variation assessment. Scientists at Roche Nimblegen have shown that the EZ DNA Methylation-Lightning™ kits work quite well in this system.

Kit Recommendation: EZ DNA Methylation-Lightning™ Kit



Get a Zymo Sample Kit. Test out any of the Zymo Research bisulfite conversion kits without risk. Just fill out the Sample Kit form, found on each product page, to request a trial version of your selected product.



EZ DNA Methylation™ kits are recommended by:

- 1 Illumina Infinium® and GoldenGate® assay systems
- 2 Sequenom MassARRAY® EpiTYPER®
- 3 Epicentre EpiGnome™ Methyl-Seq kit
- 4 Agilent Technologies SureSelect^{XT} Methyl-Seq Target Enrichment System
- 5 Roche NimbleGen SeqCap Epi Enrichment System

SUMMARY

Bisulfite conversion has been, and continues to be a cornerstone of DNA methylation analysis methodologies. Although there are some considerable challenges associated with bisulfite treatment, they can be overcome relatively easily using the techniques that we have discussed in this guide. If you have additional questions, or would like to learn more about the application of bisulfite conversion kits in your research, please contact Zymo Research Technical Support at tech@zymoresearch.com

References

1. Shen L, Guo Y, Chen X, Ahmed S, Issa JP. Optimizing annealing temperature overcomes bias in bisulfite PCR methylation analysis. *Biotechniques*. 2007 Jan;42(1):48, 50, 52 passim.

PRODUCTS

BISULFITE CONVERSION

Product	Size	Cat. No.
EZ DNA Methylation™ Kit	50 Preps.	D5001
	200 Preps.	D5002
EZ-96 DNA Methylation™ Kit (Shallow Well Format)	2x96 Preps.	D5003
EZ-96 DNA Methylation™ Kit (Deep Well Format)	2x96 Preps.	D5004
EZ-96 DNA Methylation™ MagPrep	4x96 Preps.	D5040
	8x96 Preps.	D5041
EZ DNA Methylation-Gold™ Kit	50 Preps.	D5005
	200 Preps.	D5006
EZ-96 DNA Methylation-Gold™ Kit (Shallow Well Format)	2x96 Preps.	D5007
EZ-96 DNA Methylation-Gold™ Kit (Deep Well Format)	2x96 Preps.	D5008
EZ-96 DNA Methylation-Gold™ MagPrep	4x96 Preps.	D5042
	8x96 Preps.	D5043
EZ DNA Methylation-Direct™ Kit	50 Preps.	D5020
	200 Preps.	D5021
EZ-96 DNA Methylation-Direct™ Kit (Shallow Well Format)	2x96 Preps.	D5022
EZ-96 DNA Methylation-Direct™ Kit (Deep Well Format)	2x96 Preps.	D5023
EZ-96 DNA Methylation-Direct™ MagPrep	4x96 Preps.	D5044
	8x96 Preps.	D5045
EZ DNA Methylation-Lightning™ Kit	50 Preps.	D5030
	200 Preps.	D5031
EZ-96 DNA Methylation-Lightning™ Kit (Shallow Well)	2x96 Preps.	D5032
EZ-96 DNA Methylation-Lightning™ Kit (Deep Well)	2x96 Preps.	D5033
EZ-96 DNA Methylation-Lightning™ MagPrep	4x96 Preps.	D5046
	8x96 Preps.	D5047

UNIVERSAL METHYLATED DNA STANDARDS

Product	Size	Cat. No.
Human Methylated & Non-Methylated DNA Set	1 Set	D5014
Human WGA Methylated & Non-Methylated DNA Set	1 Set	D5013
Universal Methylated Human DNA Standard	20 Rxns. (5 µg)	D5011
Universal Methylated Mouse DNA Standard	20 Rxns. (5 µg)	D5012

5-MC AND 5-HMC QUANTIFICATION

Product	Size	Cat. No.
5-mC DNA ELISA Kit	1x96 Preps.	D5325
	2x96 Preps.	D5326
Quest 5-hmC™ DNA ELISA Kit	1x96 Preps.	D5425
	2x96 Preps.	D5426

Phone: (949) 679-1190

Email: info@zymoresearch.com

www.zymoresearch.com

These products are intended for Research Use Only. Not for use in diagnostic procedures. Information, descriptions, and specifications in this publication are subject to change without notice.


ZYMO RESEARCH
The Beauty of Science is to Make Things Simple

THE
Epigenetics
COMPANY™